

MicroResp™: A VERSATILE METHODOLOGY FOR MEASURING ENVIRONMENTAL QUA

Clare M. Cameron¹, Colin D. Campbell^{1,2}, Stephen J. Chapman¹, Duncan White¹

¹Macaulay Land Use Research Institute, Aberdeen, AB15 8QH, UK; ²Swedish Agricultural Sciences University, Sweden • c.cameron@macaulay.ac.uk


INTRODUCTION: A primary role of soil is the degradation of organic matter and xenobiotic compounds that underpin nutrient cycling and environmental protection. Changes in the diversity of the soil microbial community is often used to assess soil health and one way this can be achieved is by measuring their respiration induced by adding a range of carbon substrates.

METHODOLOGY: What is MicroResp™?

- a microplate-based respiration system that overcomes the disadvantages of both classic carbon source growth tests (e.g. Biolog™) and the space and time needed for conventional respiration methods.
- measures the immediate (0-6h) substrate induced response reflecting activity rather than growth.
- can be used for long term degradation kinetics by seguential reads.

MicroResp™ measures the CO₂ evolved from soil and carbon substrate utilisation in two ways:





RESULTS:

(293 samples in 3 weeks) from a national soil sampling programme to relate respiration to land use and soil C using the colorimetric method (Fig. 1)

i) High through-put testing of large numbers of soils ii) Testing the degradation of (14C-labelled) xenobiotics in soils challenged with elevated levels of Cu and Zn (Fig.2)




Figure 1. Relationship of %C and respiration of soils of different land use




Figure 2. Degradation of xenobiotics added to sewage sludge treated soils

CONCLUSIONS: The ability to vary the C sources and the soil/sediment test conditions mean that the method is very versatile for testing a range of environmental quality issues at field and regional scales. Recent applications in testing decomposition processes in waters and sediments are also now being explored.